















Features

- · Verification class III (verification is optional)
- · Approved as a medical device according to 93/42/EEC or regulation (EU) 2017/745
- · This chair scale is the ideal measuring instrument for retirement homes, rehabilitation centres and clinics with lots of patients, and can weigh overweight people up to 250 kg
- 1 Mobile version with two steerable rollers and particularly convenient locking brakes at the rear
- · With its four wheels, this chair scale provides maximum mobility in bringing the scale to the patient. This ensures more efficient use of time for the clinic staff and greater safety for the patients who can be weighed in their familiar environment
- · The rollers, with their large diameter, make it easier to get over door thresholds, edges and across the gaps into elevators
- · For fragile patients, the extra wide, comfortable, ergonomically-optimised seat offers safe, secure seating during weighing
- · Because of the clearly contrasting black colour of the seating surface, footrests and armrests ideally suitable for dementia patients
- 2 Two foldable armrests and footrests make transfers into the chair easier. Ideal for overweight patients or barrier-free use, e.g. for transfers from a bed to the scale
- 3 Ergonomically positioned carrying handles

- · Hold function: While weighing patients that are unable to sit still, a mean average weight value is determined. This allows for sufficient time to attend to the patient, and then get a weight reading
- · BMI function to determine underweight/normal weight/surplus weight
- · High resolution readability: readability [d] can be increased by one decimal place for 5 sec. by the touch of a key
- · Protective working cover included with delivery

Technical data

- S Large backlit LCD display, digit height 25 mm
- · Dimensions of display device W×D×H 200×130×60 mm
- · Dimensions of seating surface W×D 390×360 mm
- Overall dimensions W×D×H 625×990×985 mm
- Rechargeable battery pack integrated, as standard, operating time up to 40 h without backlight, charging time approx. 12 h
- · Battery operation possible, 6×1.5 V AA not included, operating time up to 20 h
- · 4 External mains adapter as standard, with strain-relief and disconnector plug to protect the power supply components
- · Net weight approx. 24 kg

Accessories

- Dractical mains adapter pouch to store the mains adapter which is supplied as standard. Fitting option using two Velcro fasteners, KERN MCC-A01
- Protective working cover over the display device, scope of delivery: 5 items, KERN MBC-A06S05
- · Cleaning cloths, alcohol-free cloths for disinfectant wiping, quick acting, based on modern quaternary ammonium compounds and effective against papova viruses. Particularly gentle on materials, and very well suited for disinfecting products which are sensitive to alcohol. Fulfill the legal requirements for occupational safety in accordance with TRGS 525/540. Packaging contents 100 pcs., size 20×22 cm per cloth, KERN MYC-01

*Within the EU, official verification (conformity assessment according to NAWI 2014/31/EU) is mandatory by law for scales that are intended for use as a medical device. Please add this to your order. We require the location of use and the post code for the verification















| Model | Weighing range | Readout | Verification value | Mandatory by law Verification |
|--------------|----------------|-----------|--------------------|--------------------------------|
| KERN | [Max] kg | [d] kg | [e] kg | MIII) KERN |
| MCC 250K100M | 250 | 0,1 | 0,1 | 965-129 |